

INTRODUCING THE

CellFX[®] System

POWERED BY NANO-PULSE STIMULATION[™] TECHNOLOGY

A new multi-application platform for difficult-to-treat benign skin lesions. The CellFX System uses an innovative non-thermal technology that uniquely works at the cellular level where the skin problem resides. The CellFX procedure is a quick and elegant method to clear cellular-based lesions while minimizing risk of scarring.^{1,2}



The Wow is in the How

The CellFX System uses an entirely new energy modality called Nano-Pulse Stimulation (NPS[™]) technology that uses ultrafast, high-amplitude electrical energy pulses to clear cellular structures while sparing surrounding non-cellular dermal foundation.

CELL-SPECIFIC

- Affects cells and only cells
- No effect on non-cellular collagen foundation

UNIQUE MARGIN CONTROL


- Ability to reach cellular structures in the treatment area per tip selection
- Naturally seeks where abnormal cells stop and healthy non-cellular dermis begins

NON-THERMAL

- Non-thermal energy completely unlike any other modality, including lasers, cryo and RF
- Minimizes risk of thermal damage, helping to provide favorable healing response^{1,2}

INDICATED FOR COMMON BENIGN LESIONS SUCH AS:

Seborrheic Keratoses



Sebaceous Hyperplasia



Cutaneous, Non-Genital Warts



1. Hruza G, Zelickson B, Selim M, Rohrer T, Newman J, et al. Safety and Efficacy of Nanosecond Pulsed Electric Field Treatment of Seborrheic Keratoses, *Dermatologic Surgery*. September 2020 - Volume 46 - Issue 9 - p 1183-1189 doi: 10.1097/DSS.0000000000002278 2. Munavalli G, Zelickson B, Selim M, Kilmer S, Rohrer T, Newman J, et al. Safety and Efficacy of Nanosecond Pulsed Electric Field Treatment of Sebaceous Gland Hyperplasia. *Dermatologic Surgery*. 2020; 46: 803-809 3. Nuccitelli R, LaTowsky BM, Lain E, Munavalli G, Loss L, Ross EV, Jauregui L and Knape WA (2021). Safety and Efficacy of Nano-Pulse Stimulation Treatment of Non-Genital, Cutaneous Warts (Verrucae). *Lasers in Surgery and Medicine*. <https://doi.org/10.1002/lsm.23423>

CellFX System

CellFX Console

Designed for ease of use, portability and simplicity

Tunable procedure settings to address various treatment needs

Utilization-based business model enabled by its internet connectivity to a proprietary Pulse Biosciences cloud network



CellFX Handpiece

Multi-use handpiece lets you control the delivery of pulses ("cycles")



CellFX Treatment Tips

Range of single-patient use treatment tips to accommodate a variety of lesion sizes


Automated treatment settings based on treatment tip being utilized

	Area of Treatment Tip
 1.5	
 2.5	
 5	
 7.5	
 10	

CellFX Marketplace Customer Portal

One-stop destination for CellFX System business services:

- View CellFX System usage history and session summaries
- Purchase treatment tips, cycle units and replacement parts
- View purchase history and more



The CellFX® System is indicated for dermatologic procedures requiring ablation and resurfacing of the skin for the reduction, removal, and/or clearance of cellular-based benign lesions.

DISTRIBUTOR INFORMATION HERE



0044

Manufactured by
Pulse Biosciences, Inc.
3957 Point Eden Way
Hayward, CA 94545 USA



Emergo Europe B.V.
Prinsessegracht 20
2514 AP The Hague
The Netherlands

